





## How AnteoBind has helped with the Study?

RBDCoV-ACE2 is fast, suitable for high throughput assays and highly reproducible. In instances where the standard covalent conjugation strategies, such as EDC/sNHS chemistry, are susceptible to specific mutations within the RBD region, AnteoTech's Activation Kit for Multiplex Microspheres (A-LMPAKMM-400) containing AnteoBind efficiently binds mutant RBD proteins from 11 different SARS-CoV-2 variants, including the beta strain, as illustrated in Figure 2.



**Figure 2:** ACE2-RBD binding comparison of beads with different conjugation strategies. RBD proteins from wild-type, alpha and beta variants conjugated to MagPlex beads using AnteoTech's Activation Kit (a) and EDC/sNHS chemistry (b) were incubated with biotinylated ACE2 ranging from 0.1 ng/mL to 2.5 µg/mL.


During the assay, all beads are pooled to produce a bead mix, which is then incubated with patient serum diluted in assay buffer containing biotinylated ACE2. RBD-bound ACE2 is then detected using phycoerythrin-labelled streptavidin and MFI values are normalized. RBDCoV-ACE2 generates highly reproducible data with minimal variation (as determined by low intra and inter assay variation, %CVs, Figure 3).




**Figure 3:** Results of intra assay precision (a) and inter assay precision (b) of RBDCoV-ACE2 analyzing ACE2 binding inhibition to wild-type RBD (in %) of five serum samples (vaccinated (n=4, Vac1-Vac4), infected (n=1, Inf1)). Adapted from (1).

## Why use AnteoBind as an alternative technology?

From a clinical perspective, the speed and ease of incorporating new variants into the assay compared to the standard virus neutralization tests make RBDCoV-ACE2 ideal for screening how ACE2 binding inhibition changes for emerging variants. The use of AnteoBind containing Activation kit for Multiplex Microspheres illustrates good batch-to-batch reproducibility for serum ACE2 binding inhibition on the Luminex platform.



**Figure 4:** The use of AnteoBind illustrates a high degree of batch-to-batch reproducibility. Two batches (batch 1 and 2) of beads were conjugated with wild-type RBD using AnteoTech's Activation Kit for Multiplex Microspheres. Correlated are the ACE2 binding inhibitions (in %) of serum samples from COVID-19 patients (n=62) with varying disease severity for both bead batches.



## AnteoBind Offers Flexible Solutions to Time Consuming Problems

AnteoTech recognises some key protein conjugation roadblocks in the Life Science industry. We believe resources and on-time delivery of products are critical to the success of healthcare products during and beyond the pandemic.

Please contact us at [sales@anteotech.com](mailto:sales@anteotech.com) for further information on our AnteoBind products.

For further information on RBDCoV-ACE2, please contact: [Alex.Dulovic@nmi.de](mailto:Alex.Dulovic@nmi.de)

## References

1. Junker, D., Dulovic, A., Becker, M. et al. COVID-19 patient serum less potently inhibits ACE2-RBD binding for various SARS-CoV-2 RBD mutants. *Sci Rep* 12, 7168 (2022). <https://doi.org/10.1038/s41598-022-10987-2>  
Licensed under <http://creativecommons.org/licenses/by/4.0/>

**AnteoTech Ltd**  
E [sales@anteotech.com](mailto:sales@anteotech.com)  
T +61 7 3219 0085  
[www.anteotech.com](http://www.anteotech.com)