

PCT関連 取り扱い製品



抗体

サプライヤー	製品コード	製品名	タイプ	ホスト	バッファー	備考 (検証済みアプリケーション等)	ペア抗体
Meridian Life Science	MERE86813M	Monoclonal Antibody to Human Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	MERE86412M (Capture) MERE86412M (Detection)
	MERE86420B	Monoclonal Antibody to Human Procalcitonin Biotin Conjugated	Monoclonal	Mouse	20 mM Tris-HCl, 150 mM Sodium Chloride, pH 7.5	ELISA, Western Blot	MERE86561M (Capture)
	MERE01291M	Monoclonal Antibody to Human Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	MERE86494M (Detection)
	MERE86494M	Monoclonal Antibody to Human Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	MERE01291M (Capture) MERE86561M (Capture)
	MERE86420M	Monoclonal Antibody to Human Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	MERE86561M (Capture)
	MERE01308M	Monoclonal Antibody to Human Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	MERE86241M (Detection)
	MERE01342M	Monoclonal Antibody to Procalcitonin/Katacalcin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.0	Lateral Flow	MERE01340M (Detection)
	MERE01341M	Monoclonal Antibody to Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.0	Lateral Flow	MERE01340M (Detection)
	MERE01340M	Monoclonal Antibody to Procalcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.0	Lateral Flow	MERE01341M (Capture) MERE01342M (Capture)
	MERE86412M	Monoclonal Antibody to Human Procalcitonin/Katacalcin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA	MERE86813M (Detection)
	MERE86224M	Monoclonal Antibody to Human Procalcitonin/Calcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	
	MERE86131M	Monoclonal Antibody to Human Procalcitonin/Calcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	
	MERE01329M	Monoclonal Antibody to Human Calcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA	
	MERE01253M	Monoclonal Antibody to Procalcitonin/Calcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	
MERE86561M	Monoclonal Antibody to Human Procalcitonin/Calcitonin	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA, Western Blot	MERE86420M (Detection) MERE86494M (Detection)	
BiosPacific	BPA-A91062	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	PBS	-	BPA-A91064
	BPA-A91064	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	PBS	-	BPA-A91062
	BPA-A91066	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	10mM PBS, 150mM NaCl	CLIA, CMIA	BPA-A91068
	BPA-A91068	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	10mM PBS, 150mM NaCl	CLIA, CMIA	BPA-A91066
	BPA-A91076	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	PBS, pH 7.2	-	BPA-A91077
	BPA-A91077	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	-	-	BPA-A91076
	BPA-A91078	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	PBS, pH 7.2	-	BPA-A91079
	BPA-A91079	Procalcitonin (PCT) Monoclonal Antibody, Purified	Monoclonal	Mouse	PBS, pH 7.2	-	BPA-A91078
	BPA-A91093	Procalcitonin (PCT) Monoclonal Antibody	Monoclonal	Mouse	PBS, pH 7.4	-	BPA-A91094
	BPA-A91094	Procalcitonin (PCT) Monoclonal Antibody	Monoclonal	Mouse	PBS, pH 7.4	-	BPA-A91093
Aalto Bio Reagents	ABR-BF1230	Mouse Monoclonal Antibody to Human Calcitonin 0511045	Monoclonal	Mouse	PBS buffer, pH 7.4	-	-
	ABR-BF1234	Mouse Monoclonal Antibody to Human Calcitonin 0611055	Monoclonal	Mouse	PBS buffer, pH7.4	-	(capture-detection): Clone 0511045 - clone 0911065 Clone 0611055 - clone 0511045
	ABR-BF2363	Procalcitonin, Affinity Purified (Goat)	Polyclonal	Goat	Phosphate Buffered Saline, pH 7.4	-	-
EastCoast Bio	ECB-HM563	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow, CLIA	ECB-HM564 (detection)
	ECB-HM364	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, WB	ECB-HM365 (capture)
	ECB-HM243	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl	ELISA, Lateral Flow	ECB-HM244 (capture)
	ECB-HM698	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, CLIA	ECB-HM699 (detection)
	ECB-HM564	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow, CLIA	ECB-HM563 (capture)
	ECB-HM365	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, WB	ECB-HM364 (detection)
	ECB-HM244	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl	ELISA, Lateral Flow	ECB-HM243 (detection)
	ECB-HM699	Procalcitonin (PCT) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, CLIA	ECB-HM698 (capture)
	ECB-HM778	Procalcitonin (PCT-calcitonin) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, CLIA	ECB-HM779 (detection)
	ECB-HM081	Procalcitonin (PCT-calcitonin) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	Lateral Flow	ECB-HM078 (capture)
	ECB-HM779	Procalcitonin (PCT-katacalcin) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, CLIA	ECB-HM778 (capture)
	ECB-HM220	Procalcitonin (PCT-katacalcin) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	Lateral Flow	ECB-HM078 (capture)
	ECB-HM078	Procalcitonin (PCT-N terminal) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	Lateral Flow	ECB-HM081/ ECB-HM220 (detection)

抗原

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (純度等)
Meridian Life Science	MERA01454H	Procalcitonin, Recombinant (a.a. 1-116)	Recombinant	E. coli	20 mM Tris-HCl, 150 mM Sodium Chloride, pH 8.0 prior to lyophilization.	> 95% pure (sandwich immunoassay).
	MERA01367H	Human Procalcitonin, Recombinant	Recombinant	E. coli	20 mM Tris-HCl, pH 8.0 containing 250 mM Sodium Chloride and 50% Glycerol.	> 95% pure (SDS-PAGE).
BiosPacific	BPA-391011	Human Procalcitonin (PCT) Antigen, Recombinant	Recombinant	Mammalian	1x PBS, 50% glycerol, pH 7.4	純度 : >95%
	BPA-391010	Human Procalcitonin (PCT) Antigen, Recombinant	Recombinant	E. coli	20 mM Tris HCl, 50% glycerol, pH8.0	純度 : > 95%
Aalto Bio Reagents	ABR-BF3275	Procalcitonin (PCT)	Recombinant	E. coli	Lyophilised from 20mM Tris-HCl, pH8.0, 150mM	純度 : >95%
EastCoast Bio	ECB-LA368	Procalcitonin (PCT) Antigen, Recombinant >95%	Recombinant	E. coli	50 mM Tris-HCl, pH 8.0 with 150 mM NaCl and 10% glycerin	純度 : >95% by SDS-PAGE
Lee Biosolutions	LBS496-70	Procalcitonin (PCT)	Recombinant	E. coli	phosphate buffered saline, pH 7.4 with 50% glycerol (v/v)	純度 : > 90% (SDS-PAGE)