

Calprotectin関連 取り扱い製品

抗体

サプライヤー	製品コード	製品名	タイプ	ホスト	バッファー	備考 (検証済みアプリケーション等)	ペア抗体
Meridian Life Science	MERK01376M	Monoclonal Antibody to Human Calprotectin, L1 Protein	Monoclonal	Mouse	Lyophilized from Phosphate Buffered Saline, pH 7.4.	ELISA	MERK01373M (Capture)
	MERK01375M	Monoclonal Antibody to Human Calprotectin, L1 Protein	Monoclonal	Mouse	Lyophilized from Phosphate Buffered Saline, pH 7.4.	ELISA	-
	MERK01374M	Monoclonal Antibody to Human Calprotectin, L1 Protein	Monoclonal	Mouse	Lyophilized from Phosphate Buffered Saline, pH 7.4	ELISA	-
	MERK01373M	Monoclonal Antibody to Human Calprotectin, L1 Protein	Monoclonal	Mouse	Lyophilized from Phosphate Buffered Saline, pH 7.4	ELISA	MERK01376M (Detection)
BiosPacific	BPA-A29320	Calprotectin (S100A8/A9 Complex) Monoclonal Antibody, Purified	Monoclonal	Mouse	15mM phosphate, 0.15M NaCl, pH 7.2	sandwich EIA	BPA-A29310
	BPA-A29310	Calprotectin (S100A8/A9 Complex) Monoclonal Antibody, Purified	Monoclonal	Mouse	15mM phosphate, 0.15M NaCl, pH 7.2	Latex agglutination, Lateral Flow	BPA-A29320
	BPA-A29342	Calprotectin Monoclonal Antibody, Purified	Monoclonal	Mouse	0.015M potassium phosphate, 0.85% NaCl, pH 7.2	quantitative EIA	BPA-A29341, BPA-A29343, BPA-A29344
	BPA-A29344	Calprotectin Monoclonal Antibody	Monoclonal	Mouse	0.015M potassium phosphate, 0.85% NaCl, pH 7.2	quantitative EIA	BPA-A29341, BPA-A29342
	BPA-A29343	Calprotectin Monoclonal Antibody, Purified	Monoclonal	Mouse	0.015M potassium phosphate, 0.85% NaCl, pH 7.2	quantitative EIA	BPA-A29341, BPA-A29342
	BPA-A29341	Calprotectin Monoclonal Antibody	Monoclonal	Mouse	0.015M potassium phosphate, 0.85% NaCl, pH 7.2	quantitative EIA	BPA-A29342, BPA-A29343, BPA-A29344
EastCost Bio	ECB-HM319	Calprotectin Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM320 (detection)
	ECB-HM645	Calprotectin Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM644 (capture)
	ECB-HM632	Calprotectin Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.4	ELISA, Lateral Flow, TIA	ECB-HM632 (self-pairs against calprotectin tetramer)
	ECB-HM320	Calprotectin Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM319 (capture)
	ECB-HM644	Calprotectin Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM645 (detection)
Arotec Diagnostics	AROAMC04	Anti-Calprotectin Monoclonal IgG	Monoclonal	Mouse	Phosphate buffered saline	Western blot, Direct ELISA	-
	AROABC04	Anti-Calprotectin IgG	Polyclonal	Goat	Phosphate buffered saline	Western blot, ELISA	-

抗原

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (純度等)
Meridian Life Science	MERA01671H	Human Calprotectin	Native	Human Neutrophils	20 mM HEPES-Sodium Hydroxide, pH 8.0, 0.4 M Sodium Chloride, 0.1 mM Phenylmethanesulfonyl fluoride (PMSF) and 0.2 mM DTT.	> 90% pure SDS-PAGE. Proprietary purification method.
	MERR01724	Human Calprotectin Recombinant	Recombinant	E. coli	Tris Buffer with salts, pH 8	> 95% Pure. Multi-step procedure including affinity chromatography
EastCost Bio	ECB-LA399	Calprotectin Antigen, Recombinant >95%	Recombinant	E.coli	Tris Buffer, pH 8.0	純度 : >95%
Yashraj Biotechnology	YASFAL-128	Human Calprotectin	Native	Human Neutrophils	20 mM HEPES, 0.4 M NaCl, pH 8.0±0.2	純度 : ≥95%
Arotec Diagnostics	AROATC-04	Calprotectin antigen	Native	Human neutrophils	20 mM HEPES-NaOH, pH 8.0; 0.4 M NaCl; 0.1 mM PMSF; 0.2 mM DTT; 0.0075% (w/v) Bromonitrodioxane; 0.0075% (w/v) 2-Methyl-4-isothiazolin-3-one	純度 : >90%