

IL-6関連 取り扱い製品

抗体

サプライヤー	製品コード	製品名	タイプ	ホスト	バッファー	備考 (検証済みアプリケーション等)	ペア抗体
Meridian Life Science	MERH01408M	Monoclonal Antibody to Interleukin-6 (IL-6), Recombinant	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA	MERH01407M
	MERH01407M	Monoclonal Antibody to Interleukin-6 (IL-6), Recombinant	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	ELISA	MERH01408M
BiosPacific	BPA-A10301	IL-6 Monoclonal Antibody, Purified	Monoclonal	Mouse	PBS, pH 7.4	This antibody does not cross-react with rmIL-6, rhOSM, rhLIF, rhIL-11, rhgp130 or rhCNTF.	-
	BPA-A10304	IL-6 Monoclonal Antibody	Monoclonal	Mouse	0.015 M potassium phosphate buffer, 0.85% NaCl, pH 7.2	-	-
	BPA-A10305	IL-6 Monoclonal Antibody	Monoclonal	Mouse	0.015 M potassium phosphate buffer, 0.85% NaCl, pH 7.2	-	-

抗原

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (純度等)
BiosPacific	BPA-J10311	Human IL-6 Antigen	Recombinant	CHO cell	Lyophilized from PBS	-
	BPA-J10312	Human IL-6 Antigen, Recombinant	Recombinant	E. coli	10 mM sodium phosphate, 150mM NaCl, 25% glycerol, 1mM AEBSF, pH 7.4 ± 0.2	-
	BPA-J10315	Human IL-6 Protein, Recombinant	Recombinant	E. coli	PBS + 1M NaCl, pH 7.4 with 50ug BSA per 1ug as carrier protein	-
Lee Biosolutions	LBS341-36	Interleukin-6, human	Recombinant	E. coli	sodium chloride, pH 8.0 with 50% glycerol (v/v)	純度 : ≥ 90% (SDS-PAGE)
Yashraj Biotechnology	YASFRL6-87	Recombinant IL-6 (rIL-6)	Recombinant	E. coli	10 mM sodium phosphate, 150mM NaCl, 25% glycerol, 1mM AEBSF, pH 7.4 ± 0.2	純度 : ≥ 97%