

NGAL関連 取り扱い製品



抗体

サプライヤー	製品コード	製品名	タイプ	ホスト	バッファー	備考 (検証済みアプリケーション等)	ベア抗体
Meridian Life Science	MERK01392C	Chicken Antibody to Neutrophil Gelatinase-Associated Lipocalin (NGAL)/Lipocalin-2	Polyclonal	Chicken Eggs	0.1 M Sodium Phosphate, 0.15 mM Sodium Chloride, pH 7.17	-	-
BiosPacific	BPA-A05010	Lipocalin-2/NGAL Monoclonal Antibody	Monoclonal	Rat	PBS, 5% Trehalose, pH 7.4	recombinant human (rh) Lipocalin-1, rhMMP-9, rmLipocalin-2, rrLipocalin-2 との交差反応性なし	-
	BPA-A05011	Lipocalin-2/NGAL Monoclonal Antibody	Monoclonal	Mouse	PBS, pH 7.4	recombinant human (rh) Lipocalin-1, rmLipocalin-2との交差反応性なし、human Lipocalin-2/MMP-9を認識	-
	BPA-A05012	NGAL Monoclonal Antibody	Monoclonal	Mouse	0.015 M Potassium Phosphate buffer containing 0.85% NaCl, pH 7.2	-	BPA-A05013
	BPA-A05013	NGAL Monoclonal Antibody	Monoclonal	Mouse	0.015 M Potassium Phosphate buffer containing 0.85% NaCl, pH 7.2	-	BPA-A05012, BPA-A05014, BPA-A05015
	BPA-A05014	NGAL Monoclonal Antibody, Purified	Monoclonal	Mouse	0.015 M Potassium Phosphate buffer containing 0.85% NaCl, pH 7.2	human NGAL との反応性確認済み	BPA-A05013
	BPA-A05015	NGAL Monoclonal Antibody, Purified	Monoclonal	Mouse	0.015 M Potassium Phosphate buffer containing 0.85% NaCl, pH 7.2	human NGAL との反応性確認済み	BPA-A05013
	BPA-A05016	NGAL Monoclonal Antibody	Monoclonal	Mouse	PBS, pH 7.4	-	BPA-A05017
BPA-A05017	NGAL Monoclonal Antibody	Monoclonal	Mouse	PBS, pH 7.4	-	BPA-A05016	
Aalto Bio Reagents	ABR-BF1363	Mouse Monoclonal Antibody to Human Neutrophil Gelatinase-Associated Lipocalin (NGAL)	Monoclonal	Mouse	PBS buffer, pH 7.4	ヒトNGALモノマー、ホモダイマーを認識 Clone 0303011, 0303021, 0303051 は NGAL NGAL/MMP-9 ヘテロダイマーも認識	Clone 0303021 - Clone 0303011 Clone 0303021 - Clone 0303051 ABR-BF1365: Clone 0403061 - Clone 0303031 Clone 0403071 - Clone 0303011 Clone 0403071 - Clone 0303041 Clone 0403071 - Clone 0303051
	ABR-BF1365	Rabbit Monoclonal Antibody to Human Neutrophil Gelatinase-Associated Lipocalin (NGAL)	Monoclonal	Rabbit	PBS buffer, pH 7.4	ヒトNGALモノマー、ホモダイマーを認識 Clone 0403071 は NGAL NGAL/MMP-9 ヘテロダイマーも認識	ABR-BF1363: Clone 0403061 - Clone 0303031 Clone 0403071 - Clone 0303011 Clone 0403071 - Clone 0303041 Clone 0403071 - Clone 0303051
EastCoast Bio	ECB-HM306	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM307
	ECB-HM307	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM306
	ECB-HM308	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM309
	ECB-HM309	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.0	ELISA, Lateral Flow	ECB-HM308
Lee Biosolutions	LBSGNGAL-80A-Z	anti-Human NGAL	Polyclonal	Goat	phosphate buffered saline, pH 7.2	-	-

抗原

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (純度等)
Meridian Life Science	MER9278	NGAL Antigen, Recombinant	Recombinant	CHO cells	Phosphate Buffered Saline	-
BiosPacific	BPA-J05000	NGAL Antigen	Native	Human Neutrophils	50mM sodium acetate, 0.25M sodium chloride, 22% Glycerol, pH 4.5 ± 0.2	-
Aalto Bio Reagents	ABR-BF3369	Recombinant Human Human Neutrophil Gelatinase-Associated Lipocalin (NGAL)	Recombinant	mammalian cell	freeze-dried from 20mM Tris-HCl buffer, pH 7.5, containing 0.15 sodium chloride and 0.1% CHAPS	-
EastCoast Bio	ECB-LA312	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antigen, Recombinant >98%	Recombinant	E. coli	10 mM Phosphate Buffered Saline, pH 7.0	Purified Recombinant Antigen (Full length with His-tag)
Lee Biosolutions	LBS342-49	Neutrophil Gelatinase-Associated Lipocalin	Native	Human Leukocytes	sodium acetate, sodium chloride and glycerol	純度: ≥90% (SDS-PAGE)
	LBS342-49R	Neutrophil Gelatinase-Associated Lipocalin	Recombinant	<i>Pichia pastoris</i> (Yeast)	Lyophilized powder from 0.02 M ammonium bicarbonate	純度: ≥96% (SDS-PAGE)
Yashraj Biotechnology	YASFRMGAL-33	Recombinant Neutrophil Gelatinase Associated Lipocalin (rNGAL-mutant)	Recombinant	E. coli	50 mM Sodium Phosphate, 0.15M Sodium Chloride, pH7.4±0.2 with 50% Glycerol	純度: ≥95% (SDS-PAGE)
	YASFRWAL-34	Recombinant Neutrophil Gelatinase Associated Lipocalin (rNGAL-wild)	Recombinant	E. coli	50 mM Sodium Phosphate, 0.15M Sodium Chloride, pH7.4±0.2 with 50% Glycerol	純度: ≥95% (SDS-PAGE)