

FSH関連 取り扱い製品



抗体

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (検証済みアプリケーション等)	ヘア抗体
Meridian Life Science	MERMCF02-011	Monoclonal Antibody to Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF02-264 MERMAF02-275
	MERMAF02-275	Monoclonal Antibody to Human Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMCF02-011 MERMAF02-243
	MERMAF02-264	Monoclonal Antibody to Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMCF02-011 MERMAF02-243 MERMAF02-145
	MERMAF02-243	Monoclonal Antibody to Human Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF02-264 MERMAF02-275
	MERMAF02-155	Monoclonal Antibody to Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF02-264
	MERMAF02-145	Monoclonal Antibody to Human Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF02-264
	MERE86123M	Monoclonal Antibody to Follicle Stimulating Hormone (FSH-b), beta Subunit	Monoclonal	Mouse	10 mM Tris, 0.15 M Sodium Chloride, pH 7.5	ELISA, WB	MERE86051M
	MERE87701M	Monoclonal Antibody to Human Follicle Stimulating Hormone (FSH) Intact, cl Epitope	Monoclonal	Mouse	Phosphate Buffered Saline	ELISA, RIA	
MERE45709M	Monoclonal Antibody to Follicle Stimulating Hormone (FSH)	Monoclonal	Mouse	0.015 M Potassium Phosphate Buffer, 0.15 M Sodium Chloride, pH 7.2			
BiosPacific	BPA-A18030	Human FSH Monoclonal Antibody, Purified	Monoclonal	Mouse	10 mM phosphate, 150 mM NaCl, pH 7.4		
	BPA-A18052	Human FSH Monoclonal Antibody, Purified	Monoclonal	Mouse	0.015 M potassium phosphate buffer, 0.85% NaCl, pH 7.2		BPA-A18081 (monoclonal Alpha-hCG : hCGと共通の構造を認識する抗体)
	BPA-A18053	Human FSH Monoclonal Antibody, Purified	Monoclonal	Mouse	0.015 M potassium phosphate buffer, 0.85% NaCl, pH 7.2		BPA-A18081 (monoclonal Alpha-hCG : hCGと共通の構造を認識する抗体)
Aalto Bio Reagents	ABR-BF1343	Mouse Monoclonal Antibody to β -Follicle Stimulating Hormone (β -FSH)	Monoclonal	Mouse	10mM Tris, 0.15 M NaCl, pH 7.5	sandwich immunoassay	ABR-BF1348 (monoclonal to hCG : hCGと共通の構造を認識する抗体)
EastCoast Bio	ECB-HM670	Follicle Stimulating Hormone α (FSH α) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA, WB	ECB-HM669 (capture)
	ECB-HM671	Follicle Stimulating Hormone α (FSH α)- HRP Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA, WB	ECB-HM669 (capture)
	ECB-FA311	Follicle Stimulating Hormone β (FSH β) Antibody	Monoclonal	Mouse	0.05 M Phosphate Buffered Saline, pH 7.4	ELISA	
	ECB-HM669	Follicle Stimulating Hormone β (FSH β) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA, WB	ECB-HM670 or ECB-HM671 (detection)

抗原

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (純度等)
Meridian Life Science	MERR01716	Follicle Stimulating Hormone (FSH) Recombinant	Recombinant	CHO Cell Line	Phosphate Buffered Saline, pH 7.4 \pm 0.1	
	MERH6F02-323	FSH > 98% Follicle Stimulating Hormone (FSH)	Native	Human pituitary glands	Lyophilized from 50 mM Ammonium Bicarbonate	> 98%
	MERA01283H	Human Follicle Stimulating Hormone (FSH) > 98% Pure	Native	Human urine of post-menopausal women	Lyophilized from water containing no additives	> 98%
BiosPacific	BPA-J18040	Human Follicle Stimulating Hormone (FSH) Antigen	Native	Human tissue	Lyophilized from 10mM NH ₄ HCO ₃	\geq 95%
	BPA-J18050	Human Follicle Stimulating Hormone (FSH) Antigen, Recombinant	Native	Chinese hamster ovary cells	Lyophilized from a 0.2 μ filtered solution in PBS with 5% Trehalose, pH 7.2	>95% by SDS-PAGE
Aalto Bio Reagents	ABR-CW3310	Part Purified Human Follicle Stimulating Hormone (\geq 75%)	Native	Human pituitary glands	Lyophilized from 50mM ammonium bicarbonate	\geq 75%
	ABR-CW3360	Highly Purified Human Follicle Stimulating Hormone (\geq 95%)	Native	Human pituitary glands	Lyophilized from 50mM ammonium bicarbonate	\geq 95%
	ABR-CW3304	Affinity Purified Human Follicle Stimulating Hormone (\geq 95%)	Native	Human pituitary glands	Lyophilized from 50mM ammonium bicarbonate	\geq 95%
Lee Biosolutions	LBS996-11	Follicle Stimulating Hormone (FSH)	Native	Human Pituitary Glands	Lyophilized from 50 mM ammonium bicarbonate.	> 95%
	LBS996-10	Follicle Stimulating Hormone (FSH)	Native	Human Pituitary Glands	Lyophilized from 50 mM ammonium bicarbonate.	> 50%