

プロゲステロン関連 取り扱い製品

抗体

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (検証済みアプリケーション等)	ペア抗体
Meridian Life Science	MERE82321M	Monoclonal Antibody to Progesterone	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.15	ELISA	
	MERE86222M	Monoclonal Antibody to Progesterone	Monoclonal	Mouse	Phosphate Buffered Saline, pH 7.4	Competitive assay of progesterone	
	MERE86702M	Monoclonal Antibody to Progesterone	Monoclonal	Mouse	0.15 M Sodium Chloride, 10 mM Tris, pH 7.5	ELISA	
BiosPacific	BPA-A25050	Progesterone Monoclonal Antibody	Monoclonal	Mouse	PBS, pH 7.4		
	BPA-A25049	Progesterone Monoclonal Antibody, Purified	Monoclonal	Mouse	50mM NaPo4, pH 7.2±0.2, 0.02% Proclin		
Aalto Bio Reagents	ABR-BF1345	Mouse Monoclonal Antibody to Progesterone	Monoclonal	Mouse	10mM Tris, 0.15 M NaCl, pH 7.5 または PBS, pH 7 (クローンによる)	competitive assay	(クローンによる)
EastCoast Bio	ECB-HM150	Progesterone Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA	ECB-LA251 (detection)
	ECB-HM525	Progesterone Antibody	Monoclonal	Mouse	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride	ELISA, Lateral Flow	ECB-LA330 (capture)
	ECB-P01-92-11M-P	Progesterone Antibody	Monoclonal	Mouse	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride	ELISA	ECB-P91-92-11H (detection)
	ECB-HM526	Progesterone Antibody	Monoclonal	Mouse	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride	ELISA, Lateral Flow	ECB-LA330 (capture)
	ECB-HM524	Progesterone Antibody	Monoclonal	Mouse	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride	ELISA, Lateral Flow	ECB-LA330 (capture)
	ECB-HM496	Progesterone, Canine Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA	ECB-LA251 (detection)
	ECB-HM497	Progesterone, Feline Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA	ECB-LA251 (detection)
	ECB-JP063	Progesterone-11 Antibody	Polyclonal	Rabbit	neat serum	ELISA	
Bioventix	ECB-JP060	Progesterone-17-OH Antibody	Polyclonal	Rabbit	neat serum		
	BVX-PROG.3H3	Sheep Monoclonal Antibody (SMA) to: progesterone (prog.)	Monoclonal	Sheep	PBS +0.1% Sodium azide at 7.3-7.4	ELISA	

抗原

サプライヤー	製品コード	製品名	バッファー
Meridian Life Science	MERA82606P	3-CMO-Progesterone-HRP Conjugate	Tris, pH 7.6 containing protein stabilizers.
BiosPacific	BPA-V54230	Progesterone, 11-ALP Conjugate	1:1 glycerol / H2O, 0.005MMgCl
	BPA-V56050	Progesterone, 11-BSA Conjugate	0.02-M potassium phosphate, pH 7.5
	BPA-V54220	Progesterone, 11-HRP Conjugate	1:1 glycerol / 0.015M phosphate, 0.15M NaCl, pH 7.2
EastCoast Bio	ECB-LA199	Progesterone-11-hemisuccinyl-BSA Antigen	10 mM Phosphate Buffered Saline
	ECB-LA251	Progesterone-11-hemisuccinyl-HRP Antigen	0.1M PBS with 0.1% Tween-20, 0.5% solubilized casein and dye
	ECB-LA244	Progesterone-17-OH-3-CMO-BSA Antigen	Lyophilized from water
	ECB-LA330	Progesterone-3-BSA Antigen	10 mM Phosphate Buffered Saline, pH 7.2
	ECB-LA651	Progesterone-3-BSA Antigen	10 mM Phosphate Buffered Saline, pH 7.2
Lee Biosolutions	ECB-LA197	Progesterone-3-CMO:BSA Antigen	Lyophilized from water
	LBS498-20	Progesterone antigen	Solution in 15 mM phosphate, 150 mM sodium chloride, pH 7.2 with 0.1% sodium azide.