

プロラクチン関連 取り扱い製品

抗体

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (検証済みアプリケーション等)	ペア抗体
Meridian Life Science	MERMAF10-710	Monoclonal Antibody to Human Prolactin	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF10-111
	MERMAF10-245	Monoclonal Antibody to Prolactin	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 ± 0.2 containing 150 mM Sodium Chloride	ELISA	MERMAF10-216
	MERMAF10-216	Monoclonal Antibody to Prolactin	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF10-245
	MERMAF10-145	Monoclonal Antibody to Prolactin	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF10-111
	MERMAF10-111	Monoclonal Antibody to Human Prolactin	Monoclonal	Mouse	10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride	ELISA	MERMAF10-710 MERMAF10-145
	MERE20630M	Monoclonal Antibody to Human Prolactin	Monoclonal	Mouse	0.15 M Phosphate Buffered Saline, pH 7.4	ELISA	MERE20122M
MERE20122M	Monoclonal Antibody to Human Prolactin	Monoclonal	Mouse	0.15 M Phosphate Buffered Saline, pH 7.4	ELISA	MERE20630M	
BiosPacific	BPA-A44052	Prolactin Monoclonal Antibody, Purified	Monoclonal	Mouse	0.1M PBS, pH 7.4	ELISA	BPA-A44051
	BPA-A44051	Prolactin Monoclonal Antibody	Monoclonal	Mouse	※お問合せ下さい		BPA-A44052
	BPA-A44040	Human Prolactin Monoclonal Antibody	Monoclonal	Mouse	10mM Phosphate, 150 mM NaCl, pH 7.4		BPA-A44020
	BPA-A44020	Human Prolactin Monoclonal Antibody	Monoclonal	Mouse	10mM Phosphate, 150 mM NaCl, pH 7.4		BPA-A44040
BPA-A44030	Prolactin Monoclonal Antibody, Purified	Monoclonal	Mouse	10mM phosphate, 150mM NaCl, pH 7.4			
Aalto Bio Reagents	ABR-BF1342	Mouse Monoclonal Antibody to Prolactin	Monoclonal	Mouse	PBS, pH 7.5	ELISA	各クローンの組み合わせにより可能
EastCoast Bio	ECB-FA312	Prolactin (PRL) Antibody	Monoclonal	Mouse	0.05M Phosphate Buffered Saline, pH 7.4	ELISA	ECB-HM120 (capture)
	ECB-HM664	Prolactin (PRL) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA, WB	ECB-HM663 (capture)
	ECB-HM120	Prolactin (PRL) Antibody	Monoclonal	Mouse	0.05M Phosphate Buffered Saline, pH 7.4	ELISA	ECB-FA310 or ECB-FA312 (detection)
	ECB-HM663	Prolactin (PRL) Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA, WB	ECB-HM664 or ECB-HM665 (detection)
	ECB-FA310	Prolactin (PRL) Antibody	Monoclonal	Mouse	0.05M Phosphate Buffered Saline, pH 7.4	ELISA	ECB-HM120 (capture)
	ECB-HM665	Prolactin (PRL)-HRP Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline, pH 7.5	ELISA, WB	ECB-HM663 (capture)
Lee Biosolutions	LBS996-42	Prolactin (PRL), anti-Human	Monoclonal	Mouse	Solution in phosphate buffered saline, pH 7.4 with 0.1% sodium azide.		Detection Ab: LBS996-43
	LBS996-43	Prolactin (PRL), anti-Human	Monoclonal	Mouse	Solution in phosphate buffered saline, pH 7.4 with 0.1% sodium azide.		Capture Ab: LBS996-42

抗原

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (純度等)
Meridian Life Science	MERA01467H	Human Prolactin Recombinant	Recombinant	<i>E. coli</i>	0.45% (w/v) Sodium Chloride, 0.4% (w/v) Sodium Phosphate, 3% (w/v) Trehalose, 0.3% (w/v) Arginine, 0.01% (v/v) Tween 20, pH 6.99	
	MERA01274H	Recombinant Human Prolactin (PRL)	Recombinant	<i>Pichia pastoris</i>	Lyophilized from 50 mM Ammonium Bicarbonate and 1 mM EDTA	
BiosPacific	BPA-J15120	Prolactin Antigen, Standard grade	Native	Human pituitary glands	Lyophilized powder; lyophilized from 50 mM ammonium bicarbonate	
	BPA-J15150	Prolactin Antigen, High Pure	Native	Human	Lyophilized	
	BPA-J15130	Prolactin, Recombinant Human Antigen	Recombinant	<i>E. coli</i>	Lyophilized powder; lyophilized from a 1 mg/ml solution in 10 mM sodium phosphate buffer, 50 mM NaCl, pH 8.0	
	BPA-J15140	Human Prolactin Antigen, Recombinant	Recombinant	<i>E. coli</i>	20 mM Tris-HCl, 10 mM EDTA, pH 8	
Aalto Bio Reagents	ABR-CT3341	Recombinant Prolactin	Recombinant	<i>E. coli</i>	0.45% (w/v) NaCl, 0.4% (w/v) Sodium Phosphate, 3% (w/v) Trehalose, 0.3% (w/v) Arginine, 0.01% Tween 20, pH 6.99 ±0.2	
	ABR-CW3311	Purified Human Prolactin ≥ 50%	Native	Human pituitary glands	Lyophilized from 50mM ammonium bicarbonate	
	ABR-CW3307	Purified Human Prolactin ≥ 95%	Native	Human pituitary glands	Lyophilized from 50mM ammonium bicarbonate	
Lee Biosolutions	LBS996-41	Prolactin (PRL)	Native	Human Pituitary Glands	Lyophilized from 50 mM ammonium bicarbonate.	> 98% (SDS-PAGE)
	LBS996-40	Prolactin (PRL)	Native	Human Pituitary Glands	Lyophilized from 50 mM ammonium bicarbonate.	> 50% (SDS-PAGE)