

テストステロン関連 取り扱い製品

抗体

サプライヤー	製品コード	製品名	タイプ	由来	バッファー	備考 (検証済みアプリケーション等)	ペア抗体
Meridian Life Science	MERE86902M	Monoclonal Antibody to Testosterone	Monoclonal	Mouse	10 mM Tris, 0.15 M NaCl, pH 7.5	ELISA	
	MERE45703M	Monoclonal Antibody to Testosterone	Monoclonal	Mouse	0.015 M Potassium Phosphate Buffer, 0.15 M Sodium Chloride, pH 7.2		
	MERE30412M	Monoclonal Antibody to Testosterone	Monoclonal	Mouse	Phosphate Buffered Saline	ELISA	
	MERE82721M	Monoclonal Antibody to Testosterone	Monoclonal	Mouse	0.05 M Potassium Phosphate, pH 7.15 containing 0.85% NaCl		
BiosPacific	BPA-A53070	Testosterone Monoclonal Antibody	Monoclonal	Mouse	0.1M Tris, 0.2M sodium chloride, 1mM EDTA, 10 mM sodium azide, pH 8.4		
Aalto Bio Reagents	ABR-CO2339	Rabbit anti-Testosterone	Polyclonal	Rabbit	neat serum		
EastCoast Bio	ECB-HM132	Testosterone Antibody	Monoclonal	Mouse	10 mM Phosphate Buffered Saline	ELISA	ECB-LA322 (detection)
	ECB-HM523	Testosterone Antibody	Monoclonal	Mouse	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride	ELISA, Lateral Flow	ECB-LA329 (capture)
	ECB-HM276	Testosterone Antibody	Monoclonal	Mouse	0.01M Sodium Chloride	ELISA, Lateral Flow	
	ECB-HM521	Testosterone Antibody	Monoclonal	Rat	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride	ELISA, Lateral Flow	ECB-LA329 (capture)
	ECB-HM522	Testosterone Antibody	Monoclonal	Mouse	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride	ELISA, Lateral Flow	ECB-LA329 (capture)
	ECB-JP066	Testosterone-19 Antibody	Polyclonal	Rabbit	neat serum		
Bioventix	BVX-TEST3.3C11	Sheep Monoclonal Antibodies (SMAs) to testosterone (testo)	Monoclonal	Sheep	PBS +0.1% Sodium azide at 7.2-7.4	ELISA	
	BVX-TEST3.36A3	Sheep Monoclonal Antibodies (SMAs) to testosterone (testo)	Monoclonal	Sheep	PBS +0.1% Sodium azide at 7.2-7.4	ELISA	

抗原

サプライヤー	製品コード	製品名	バッファー
Meridian Life Science	MERA01312B	Centrifuge before opening to ensure complete recovery of vial contents.	10 mM Phosphate, pH 7.0
BiosPacific	BPA-V56060	Testosterone, 3-BSA Conjugate	Lyophilized
	BPA-V56061	Testosterone, 3-HRP Conjugate	Tris with protein stabilizers, pH 7.6
Aalto Bio Reagents	ABR-CO2349	Rabbit anti-17- α -Hydroxyprogesterone	0.1M BBS buffer
	ECB-LA329	Testosterone-3-BSA Antigen	10 mM Phosphate Buffered Saline, pH 7.2
EastCoast Bio	ECB-FV505	Testosterone-3-BSA Antigen	10 mM Phosphate, pH 7.0
	ECB-LA322	Testosterone-3-HRP Antigen	0.1M PBS with 0.1% Tween-20, 0.5% solubilized casein and dye
Lee Biosolutions	LB5528-03	Testosterone-3	Solution in phosphate buffered saline, pH 7.2 with 0.1% sodium azide.
	LB5528-01	Testosterone-19	Solution in phosphate buffered saline, pH 7.2 with 0.1% sodium azide.