

PRODUCT SPECIFICATION

catalog number: J61000

product description: Alpha Feto Protein (AFP), High Purity

lot number: [lot specific]

source: Human cord blood

form: Liquid, aseptically filtered

purity: ≥95%

buffer: 50mM Tris,150mM NaCl, pH 7.2 ± 0.2

preservative: 0.09% NaN₃

molecular weight: ~66 kDa

AFP content: [lot specific] mg/mL (by Roche-cobas e411)

bioburden: ≤ 10 CFU/mL

total protein: [lot specific] mg/mL (OD₂₈₀)

manufacture date: [lot specific]

expiration date: [lot specific]

storage: 2-8°C

hazard/biohazard: This material contains sodium azide as a preservative; take

appropriate precaution.

Donor material has been tested and found negative for HBsAg, anti-HCV and anti-HIV 1&2 by FDA approved methods. Donors tested negative for viral DNA/RNA for HBV, HIV 1 and HCV by PCR. However, this material should still be handled with caution

as potentially hazardous.