



Normal Canine Serum

PRODUCT CODE : AS 2410

SHELF-LIFE : 1 year from date of receipt

SOURCE : Human serum from healthy canine (dog)

PRESENTATION : Supplied lyophilised from 10 mM sodium phosphate buffer, pH 7.2, containing 0.15 M NaCl.

CONCENTRATION : 60 mg/ml by Bradford (Biorad) assay using BSA as standard.

RECONSTITUTION : Reconstitute with indicated volume of de-ionised water. After reconstitution, add a preservative and store at +2 to +8°C, or aliquot and store frozen at $\leq -20^{\circ}\text{C}$

PACKAGING : 2ml and 10ml vials

PRESERVATIVE : None

STORAGE : +2 to + 8°C

