



Recombinant Follicle Stimulating Hormone (FSH)

PRODUCT CODE : AU 3306

SHELF LIFE : 2 years from date of manufacture when stored at -20°C

SOURCE : Recombinant antigen expressed in CHO cell culture.

PRESENTATION : The product is supplied liquid, 0.2µM filtered in Phosphate Buffered Saline, pH 7.4 +/- 0.1

PRESERVATIVE : None

FSH ACTIVITY : ≥ 2,500 IU/ml by Centaur FSH assay.

TOTAL PROTEIN : ≤ 35 OD units at 280nm

STORAGE TEMPERATURE : -20°C

SHIPPING TEMPERATURE : Product may be stored short term at +2 to +8°C for shipment.