



## Recombinant Luteinizing Hormone (LH)

**PRODUCT CODE :** AU 3335

**SHELF LIFE :** 3 years from date of manufacture

**SOURCE :** Recombinant antigen expressed in mammalian cell line

**PRESENTATION :** The product is supplied liquid, 0.2 $\mu$ m filtered in Phosphate Buffered Saline, pH 7.4 +/- 0.1

**PRESERVATIVE :** None

**LH ACTIVITY :** 750 - 950 IU/ml by Cobas LH assay.

**TOTAL PROTEIN :** <100 OD units at 280nm

**STORAGE TEMPERATURE :** -20°C

