



**Purified and Stabilised Recombinant
Streptolysin O (SLO) Antigen**

PRODUCT CODE : AW 6093-b

SHELF LIFE : 45 months from date of manufacture

SOURCE : Transformed E. Coli

STORAGE BUFFER : Tris-HCl buffered saline, pH 8.2

PURIFICATION : Metal affinity chromatography

PROTEIN CONCENTRATION : 1-2 mg/ml (by Lowry method)

TOTAL ANTIGEN ACTIVITY : >2500 IU/ml by internal indirect method

SPECIFIC ANTIGEN ACTIVITY : > 1250 IU/mg protein by internal indirect method

PRESERVATIVE : < 0.1% sodium azide

STORAGE : Frozen (-70°C)
Avoid repeated freezing and thawing