



Mouse Monoclonal Antibody to Helicobacter Pylori Flagellin

PRODUCT CODE : AZ 1256

CLONE NO. : 2510086

SHELF LIFE : 3 years from date of manufacture

IMMUNOGLOBULIN CLASS : IgG

CONCENTRATION : >2mg/ml by OD at 290nm

PRESENTATION : Purified by protein G chromatography and presented liquid in PBS, pH 7.5

IMMUNOGEN : Purified flagellin of native H. pylori

APPLICATION : Potentially suitable for ELISA, lateral Flow, ChLIA and other applications for detection of 54kDa flagellin of H. pylori. This antibody will self-pair and can be used as both capture and label elements in immunoassay.

PRESERVATIVE : 0.05% sodium azide

STORAGE : +2°C to +8°C