



## **Mouse Monoclonal Antibody to Human C-Reactive Protein**

**PRODUCT CODE :** AZ 1280

**CLONE NOS. :** 2002018  
2102028

**SHELF LIFE :** 3 years from date of manufacture

**IMMUNOGLOBULIN CLASS :** IgG

**CONCENTRATION :** >2mg/ml by OD at 280nm (extinction  
coefficient 1.4)

**PRESENTATION :** Purified by protein A chromatography and  
presented liquid in PBS, pH 7.4

**IMMUNOGEN :** Purified human CRP

**APPLICATION :** Clone 2102028 is recommend as capture and  
clone 2002018 as detection antibody in  
immunoassay development.

**PRESERVATIVE :** 0.05% sodium azide (clone 2102028)  
0.1% sodium azide (clone 2002018)

**STORAGE :** +2°C to +8°C

