



Mouse Monoclonal Antibody to Human D-Dimer

PRODUCT CODE : BF 1151

Clone Number	Immunoglobulin Subclass
2803011	IgG 1
2903021	IgG2a
3003031	IgG2a
3103041	IgG2a
1108031	IgG2b

PRESENTATION : Protein A purified immunoglobulin in PBS buffer, pH 7.4

SPECIFICITY : No cross-reaction with fibrinogen and D-Monomer for all clones. Clone 2803011 also recognises a cross-linked region of D-Dimer.

APPLICATION : All five clones are potentially suitable for coating to latex particles for latex enhanced immunoturbidimetric application.

PRESERVATIVE : 0.1% sodium azide

STORAGE : +2°C to +8°C



Church Lane,
Rathfarnham Village,
Dublin 14,
Ireland.

Tel: +353-1-4900685
Fax: +353-1-4900122

E-mail: info@aaltobioreagents.com
Web: www.aaltobioreagents.com