



Partially Purified CA 15-3

PRODUCT CODE : BF 3315

SHELF-LIFE : 12 months from date of manufacture

SOURCE : Human fluids

PRESENTATION : The product is supplied as a frozen liquid in 1M NaCl

PRESERVATIVE : 0.09% NaN₃

CA 15-3 ACTIVITY : KU/ml by EIA

CROSS REACTIVITY : Produce contains trace amounts of CA 19-9,
CA 125 and CEA.

**ANTI-HIV-1&2, ANTI-HCV
AND HBsAg :** Donor blood tested and confirmed negative.

STORAGE TEMPERATURE : -20°C
Avoid repeat freeze/thaw cycles. It is recommended to
aliquot into working volumes after first thawing.