



Mycoplasma Pneumoniae Antigen

PRODUCT CODE : BM 6062

SHELF LIFE : 10 years from date of manufacture

SOURCE : NCTC Strain FH

PRESENTATION : Enriched adhesin P1 antigen presented in phosphate solubilising buffer containing 2% n-octylglucopyranoside.

PROTEIN CONCENTRATION : mg/ml by Bradford (Biorad) assay using BSA as standard.

INACTIVATION : By detergent treatment. Confirmed by attempted growth under original culture conditions.

ANALYSIS : By SDS PAGE

PRESERVATIVE : None

STORAGE : -70°C. Sonicate before use.
Repeated freezing and thawing should be avoided.