



Recombinant Rubella Spike Ectodomain (E1-E2) Antigen

PRODUCT CODE : BM 6428

SHELF LIFE : 10 years from date of manufacture

SOURCE : Rubella virus strain HPV-77 spike glycoprotein ectodomain (E1 and E2) expressed in insect cells.

FORM : Presented as purified antigen in a proprietary buffer formulation.

PURIFICATION : by chromatography.

PROTEIN CONCENTRATION : 0.5 - 2mg/ml by absorption at 280nm (Extinction Coefficient 1.99).

ANALYSIS : SDS-PAGE.

MOLECULAR WEIGHT : 78kDa

STORAGE : -70°C.
Avoid repeated freezing and thawing.

APPLICATION : Suitable for ELISA, bead based assays and similar solid phase formats.