



Recombinant VZV Envelope Glycoprotein E Antigen

PRODUCT CODE : BM 6429

SHELF LIFE : Currently under evaluation (results due April 2021)

SOURCE : Envelope glycoprotein E ectodomain of Varicella Zoster Virus (VZV) expressed in insect cells.

FORM : Presented as purified antigen in a proprietary buffer formulation.

PURIFICATION : by chromatography.

PROTEIN CONCENTRATION : mg/ml by absorption at 280nm (protein concentration range TBC)

ANALYSIS : In process, SDS-PAGE and QC patient panel

MOLECULAR WEIGHT : 61kDa

STORAGE : -70°C.
Avoid repeated freezing and thawing.

APPLICATION : Potentially suitable for ELISA, bead based assays and similar solid phase formats.