



Recombinant SARS-CoV-2 Nucleocapsid Protein

PRODUCT CODE : CK 6404

SHELF-LIFE : 4 years from date of manufacture

SOURCE : Full-length nucleocapsid protein of SARS-CoV-2 (COVID-19) expressed in *E. coli*. Contains 6x His tag located at the C terminal.

PRESENTATION : Frozen liquid in PBS and 25mM K₂CO₃

PRESERVATIVE : 0.025% sodium azide

PURITY : 95% by SDS-PAGE. Distinct band observed at 48kDa.

PROTEIN CONCENTRATION : mg/ml (Lowry method)

STORAGE : ≤-70 °C. Avoid repeated freeze-thaw.

APPLICATON : Potentially suitable for ELISA or lateral flow type assay.

REV 02 17/07/20



E-mail: info@aaltobioreagents.com
Web: www.aaltobioreagents.com