



Rabbit monoclonal antibody to *Streptococcus pneumoniae*

PRODUCT CODE :	CM 1236
CLONE NO. :	1810286
SHELF LIFE :	2 years from date of manufacture
IMMUNOGEN :	<i>Streptococcus pneumoniae</i> cell wall polysaccharide (urinary antigen).
PRESENTATION :	Protein A purified from culture medium and presented as a liquid in 0.01M PBS, pH 7.2 containing 0.1% of sodium azide as preservative.
IMMUNOGLOBULIN SUBCLASS :	IgG1
SPECIFICITY :	Specific to cell wall polysaccharide (CWPS) of <i>Streptococcus pneumoniae</i> (urinary antigen)
PURITY :	>90% by immunoelectrophoresis using agarose
CONCENTRATION :	mg/ml by OD280nm
APPLICATION :	Use as capture antibody in ELISA, lateral flow and IFA with code CM 1247, clone 1710276, as label antibody.
STORAGE :	+2 to 8°C

