



## **Rabbit monoclonal antibody to Streptococcus pneumoniae**

**PRODUCT CODE :** CM 1247

**CLONE NO. :** 1710276

**SHELF LIFE :** 2 years from date of manufacture

**IMMUNOGEN :** Streptococcus pneumoniae cell wall polysaccharide (urinary antigen).

**PRESENTATION :** Protein A purified from culture medium and presented as a liquid in 0.01M PBS, pH 7.2 containing 0.1% of sodium azide as preservative.

**IMMUNOGLOBULIN SUBCLASS :** IgG1

**SPECIFICITY :** Specific to cell wall polysaccharide (CWPS) of Streptococcus pneumoniae (urinary antigen)

**PURITY :** >90% by immunoelectrophoresis using agarose

**CONCENTRATION :** mg/ml by OD280nm

**APPLICATION :** Use as label antibody in ELISA, lateral flow and IFA with code CM 1236, clone 1810286, as capture antibody.

**STORAGE :** +2 to 8°C



Church Lane  
Rathfarnham Village  
Dublin 14  
Ireland

Tel: +353-1-4900685  
Fax +353-1-4900122