



Mouse Monoclonal Antibody Canine C-reactive Protein

PRODUCT CODE : CZ 9012

Clone Number	Immunoglobulin Subclass
0204020	IgG1
0304030	IgG1

SHELF-LIFE : 1 year from date of manufacture

PRESENTATION : Purified by chromatography by Protein G in PBS buffer, pH 7.4

IMMUNOGEN : Recombinant Canine C-reactive protein

CONCENTRATION : A minimum of 1mg/ml, OD at 280nm

PRESERVATIVE : 0.1% sodium azide

ANALYSIS : By SDS PAGE

APPLICATION : The antibodies are potentially suitable for detection of canine CRP by lateral flow, ELISA and TIA applications .

STORAGE : -20°C
Product may be stored short term (up to one month) at +2°C to +8°C for shipping.

