



Mouse Monoclonal Antibody to Feline Serum Amyloid A

PRODUCT CODE : CZ 9069

Clone Number	Immunoglobulin Subclass
0604040	IgG1
0704050	IgG1

SHELF-LIFE : 1 year from date shipment

PRESENTATION : Purified by chromatography by Protein G in PBS buffer, pH 7.4

IMMUNOGEN : Recombinant Feline serum amyloid A

CONCENTRATION : A minimum of 1 mg/ml, OD at 280nm

PRESERVATIVE : 0.1% sodium azide

ANALYSIS : By SDS PAGE

APPLICATION : The antibodies are potentially suitable for detecting feline Serum Amyloid A by lateral flow and ELISA application.

STORAGE : $\leq -20^{\circ}\text{C}$.
Avoid repeated freezing and thawing

