



## PRODUCT SPECIFICATION

catalog number:	<b>J29400</b>
product description:	Human Calprotectin, Native
lot number:	[lot specific]
source:	Human Neutrophils
form:	Liquid
purification method:	Aseptically filtered
purity:	≥95% by SDS PAGE
buffer:	20mM HEPES, 0.4M NaCl, pH 8.0 ±0.2
activity:	2.0 ±0.4mg/ml
bioburden:	≤10 CFU/ml
total protein:	1.0 ±0.2mg/ml
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	Upon receipt: -20°C (in working aliquots to avoid freeze thaw cycles)
hazard/biohazard:	There is no known biohazard associated with this product.  Donor material has been tested and found negative for HBsAg, HIV 1&2 and, HCV by FDA approved methods and was also found negative of viral DNA/RNA for HBV, HIV1, and HCV by PCR. However, this material should still be handled with caution as potentially hazardous.