

PRODUCT SPECIFICATION

catalog number: J46290

product description: Human Epididymis Protein 4 (HE4) Antigen

lot number: [lot specific]

source: Human fluids

form: Liquid, aseptically filtered

purity: 90% SDS-PAGE

buffer: 50mM Sodium phosphate, pH 7.4 ±0.2 with 20% Glycerol

activity: [lot specific] nmoles/L

bioburden: <10 CFU/mL

contaminants:

• CA 19-9 [lot specific] • CA 125 [lot specific] IU/ml

IU/ml

CA 72-4 [lot specific]15-3 [lot specific] IU/ml

IU/ml

CEA [lot specific]

ng/ml

manufacture date: [lot specific]

expiration date: [lot specific]

storage: -20°C

hazard/biohazard: There is no known hazard or biohazard associated with this

product.

Donor material has been tested and found negative for HBsAg, HIV Ag, anti-HCV and anti-HIV 1&2 by FDA approved methods. Donors tested negative for viral DNA/RNA for HBV, HIV 1 and HCV by PCR. However, this material should still be handled with

caution as potentially hazardous.