



PRODUCT SPECIFICATION

catalog number:	J80110
product description:	Neuron Specific Enolase (NSE) Antigen, Recombinant, High Pure
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid, aseptically filtered
purification method:	Metal affinity chromatography
purity:	>95% by SDS-PAGE (major band at 50.8kDa)
buffer:	20mM sodium phosphate, 150mM NaCl, 50% glycerol, pH 7.4±0.2
preservative:	none
concentration:	[lot specific] mg/mL
bioburden:	≤10 CFU/mL
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	-20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.