

## PRODUCT SPECIFICATION

catalog number: J10740

product description: Human NT-ProBNP Antigen, Recombinant, High Purity

lot number: [lot specific]

source: E. coli

form: Liquid, aseptically filtered

purification method: Metal affinity chromatography

purity: ≥95% (by SDS-PAGE)

buffer: 20mM Tris, 150mM NaCl, 0.5M urea, 5% glycerol, 1mM PMSF,

pH 8.0 ±0.2

molecular weight: ~10.7 kDa

concentration; [lot specific] mg/mL (Roche ECLIA)

bioburden: <10 CFU/mL

notes: This protein was expressed using 76aa sequence and contained

an N-terminal His tag.

manufacture date: [lot specific]

expiration date: [lot specific]

storage: Shipped at 2-8°C. Upon arrival store at -70°C avoiding multiple

freeze thaw cycles.

hazard/biohazard: There is no known hazard or biohazard associated with this

product.