



抗体

| 製品コード | 製品名 | 免疫抗原 | タイプ | ホスト | バッファー | Matched pair | 備考 |
|------------|---|----------------------|------------|-------|--|---|--------------------|
| BPA-A10741 | NT-proBNP (12-16) Monoclonal Antibody, Purified | Native human proBNP | Monoclonal | Mouse | 10mM PBS, pH 7.4 | BPA-A10751, BPA-A10752 (as the capture) | NTproBNP(12-16)を認識 |
| BPA-A10742 | NT-proBNP (27-31) Monoclonal Antibody, Purified | Native human proBNP | Monoclonal | Mouse | 10mM PBS, pH 7.4 | BPA-A10751, BPA-A10752 (as the capture) | NTproBNP(27-31)を認識 |
| BPA-A10743 | NT-proBNP (13-27) Monoclonal Antibody, Purified | Recombinant NTproBNP | Monoclonal | Mouse | PBS, pH 7.4 | - | NTproBNP(13-27)を認識 |
| BPA-A10751 | NT-proBNP (38-44) Monoclonal Antibody, Purified | Native human proBNP | Monoclonal | Mouse | 10mM PBS, pH 7.4 | BPA-A10741, BPA-A10742 (as the conjugate) | NTproBNP(38-44)を認識 |
| BPA-A10752 | NT-proBNP (64-67) Monoclonal Antibody, Purified | Native human proBNP | Monoclonal | Mouse | 10mM PBS, pH 7.4 | BPA-A10741, BPA-A10742 (as the conjugate) | NTproBNP(64-67)を認識 |
| BPA-A10753 | NT-proBNP (39-76) Monoclonal Antibody, Purified | Recombinant NTproBNP | Monoclonal | Mouse | 1XPBS | - | NTproBNP(39-76)を認識 |
| BPA-A10754 | NT-proBNP (42-46) Monoclonal Antibody, Purified | Recombinant NTproBNP | Monoclonal | Mouse | 10 mM PBS, pH 7.4 | - | NTproBNP(42-46)を認識 |
| BPA-A10771 | NT-proBNP Monoclonal Antibody | Recombinant NTproBNP | Monoclonal | Mouse | PBS, pH 7.4 | BPA-A10772 (as the detection). | - |
| BPA-A10772 | NT-proBNP Monoclonal Antibody | Recombinant NTproBNP | Monoclonal | Mouse | PBS, pH 7.4 | BPA-A10771 (as the capture) | - |
| BPA-G10800 | NT-proBNP Polyclonal Antibody, Purified | Recombinant NTproBNP | Polyclonal | Goat | 10 mM sodium phosphate, 0.8% NaCl, pH 7.2; filtered in 0.2 µm solution | BPA-A10760, BPA-A10770 (as the capture) | NTproBNP(1-24)を認識 |

抗原

| 製品コード | 製品名 | 由来 | バッファー |
|------------|--|----------------|--|
| BPA-J10700 | Recombinant Human Pro-Brain natriuretic peptide (rNT-ProBNP) | <i>E. coli</i> | 50 mM potassium phosphate buffer, 250mM |
| BPA-J10710 | Recombinant Human Pro-Brain natriuretic peptide (rProBNP) | <i>E. coli</i> | 20 mM Tris HCl, 500 mM NaCl, pH 8.0 with |
| BPA-J10730 | Human ProBNP, Recombinant Antigen, High Purity | <i>E. coli</i> | 20mM sodium phosphate, 150mM NaCl, 0.5M urea, |
| BPA-J10740 | Human NT-ProBNP, Recombinant Antigen, High Purity | <i>E. coli</i> | 20mM Tris, 150mM NaCl, 0.5M urea, 5% glycerol, |