



PRODUCT SPECIFICATION

catalog number:	J47510
product description:	Streptolysin O (SLO), Recombinant Antigen
lot number:	[lot specific]
source:	<i>E. coli</i>
purity:	~98% by SDS Page (Commassie staining)
buffer:	20mM phosphate buffer, 150mM NaCl, 50% glycerol
SLO content:	[lot specific] HU/mg (by RBC Hemolysis)
bioburden:	≤10 CFU/mL
total protein:	[lot specific] mg/mL (OD280)
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	-20°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.