

抗体

製品コード	製品名	免疫抗原	タイプ	ホスト	バッファー	Matched pair	備考
BPA-A91140	Vitamin D Monoclonal Antibody, Purified	Vitamin-D derivative-BSA	Monoclonal	Mouse (腹水)	10 mM PBS, 150 mM NaCl, pH 7.4	-	-
BPA-A91156 BPA-A91157	Vitamin D Total (D2 & D3) Monoclonal Antibody	Vitamin-D derivative-BSA	Monoclonal	Mouse	0.1M PBS, pH 7.4	-	natural Vitamin D2 and Vitamin D3を認識
BPA-A91158	Vitamin D Total (D2 & D3) Monoclonal Antibody	Vitamin-D derivative-BSA	Monoclonal	Mouse	0.1M BBS, pH 8.0	-	natural Vitamin D2 and Vitamin D3を認識
BPA-A91160	25-OH Vitamin D3 Monoclonal Antibody, Purified	25-OH Vitamin D3-BSA	Monoclonal	Mouse	PBS, pH 7.4	-	-
BPA-A91170	Vitamin D Total (D2 & D3) Monoclonal Antibody	-	Monoclonal	Mouse	PBS, pH 7.4	-	-

抗原

製品コード	製品名	バッファー	備考
BPA-J91200	GC Binding Protein (Vitamin D Binding Protein) ヒト血漿由来	Lyophilized from a 20mM sodium phosphate, pH 7.4 with 150mM NaCl	-
BPA-V91140	Vitamin D Derivative - BSA	NaCl solution	抗体A91156, A91157, A91158と結合
BPA-V91141	Vitamin D Derivative - Biotin	Reconstitute in Ethanol	抗体A91156, A91157, A91158と結合
BPA-V91142	Vitamin D Derivative - Carboxylic Acid	Reconstitute in Ethanol	抗体A91156, A91157, A91158と結合
BPA-V91143	Vitamin D Derivative - Amino	Reconstitute in Ethanol	抗体A91156, A91157, A91158と結合
BPA-V91144	Vitamin D3 - BSA	PBS pH 7.4	-
BPA-V91145	Vitamin D3 - HRP	PBS pH 7.4, 50% glycerol, 1% BSA	-
BPA-V91145	Vitamin D3 - HRP	PBS pH 7.4, 50% glycerol, 1% BSA	-