

Product Datasheet

Anti-Human HBP Monoclonal Antibody, Azide Free Clone B-R37 855.500.000

Product Name	Anti-Human HBP Monoclonal Antibody, Azide Free Clone B-R37
Catalog Number	855.500.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Source	Recombinant human Azurocidin/Cap37/HBP
Specificity	Recognises both natural and recombinant human HBP
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells
Immunisation	Recombinant human Azurocidin/Cap37/HBP