



PRODUCT SPECIFICATION

catalog number:	A61300
product description:	IGF-1 Monoclonal Antibody
lot number:	[lot specific]
immunogen:	IGF-1
form:	Liquid
purification method:	Protein A affinity chromatography
purity:	≥95%
buffer:	PBS
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
application/notes:	RIA and EIA
manufacture date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.