



## Partially Purified Helicobacter Pylori Antigen

**PRODUCT CODE :** BM 6079

**SHELF LIFE :** 10 years from date of manufacture

**SOURCE :** Helicobacter pylori from ATCC 43504 is cultured in a solid medium. Washed harvested organisms are solubilised and semi-purified by detergent extraction and centrifugation.

**PRESENTATION :** In a 1% non-ionic detergent solution.

**INACTIVATION :** Detergent solubilisation. Infectivity confirmed negative by absence of growth under original cell culture conditions.

**TOTAL PROTEIN :** mg/ml by Bradford (Biorad) assay using bovine serum albumin as standard.

**PRESERVATIVE :** None.

**STORAGE :** -70°C. Sonicate before use.  
Repeated freeze/thawing of the product should be avoided as this will damage the immunological activity of the antigen.

