



Parvovirus B19 Antigen, Virus-Like Particles (VP2)

PRODUCT CODE : BM 6202

SHELF LIFE : 10 years from date of manufacture

SOURCE : Baculovirus infected sf9 cells containing virus-like particles from VP2 structural proteins of Parvovirus B19

PRESENTATION : PBS, pH 7.4 and 10% glycerol

PURIFICATION : By chromatography

PROTEIN CONCENTRATION : mg/ml by Bradford assay

ANALYSIS : SDS-PAGE and ELISA

STORAGE : -70°C
Repeated freezing and thawing should be avoided.

