



Recombinant p17 Treponema Pallidum

PRODUCT CODE : CK 6445

SHELF-LIFE : 3 years from date of manufacture

SOURCE : The full length p17 membrane protein of Treponema pallidum is expressed in E. coli. Contains 6x His tag located at the C terminal.

PRESENTATION : Presented as a frozen liquid in PBS with 25mM K₂CO₃.

PRESERVATIVE : 0.025% sodium azide

PURITY : >95% by SDS-PAGE. Distinct band visible at 22kDa.

MOLECULAR WEIGHT : 22kDa

PROTEIN CONCENTRATION : mg/ml (Lowry method)

STORAGE : ≤-20°C
Avoid repeated freeze-thaw.

APPLICATON : Well characterised for lateral flow type assay and potentially suitable for ELISA and other solid phase immunoassay formats. Centrifuge before opening to ensure complete recovery of vial contents



Unit A6
Nutgrove Office Park
Dublin D14 DX26
Ireland

T: + 353-1-4900 685