



Recombinant Feline Immunodeficiency Virus (FIV) p24 Antigen

PRODUCT CODE : CZ 9028-P

SHELF-LIFE : 1 year from date of shipment

SOURCE : Full length FIV p24 (a.a. 1 - 199) expressed in E. coli and containing GST tag.

FORM : Presented as purified antigen in PBS, pH 7.4

PURIFICATION : By affinity chromatography.

PROTEIN CONCENTRATION : mg/ml by OD at 280nm (extinction coefficient 1.0)

PRESERVATIVE : 0.1% sodium azide

ANALYSIS : SDS-PAGE. Distinct band visible at ~67kDa.

Immunoreactivity checked by ELISA against characterised positive sera.

STORAGE : $\leq -20^{\circ}\text{C}$. Ideally stored at -70°C .
Avoid repeated freezing and thawing.

APPLICATION : Potentially suitable for lateral flow and ELISA application.

