



Mouse Monoclonal Antibody to Feline Leukemia Virus

PRODUCT CODE : CZ 9111

Clone Number	Immunoglobulin Subclass
1803080	IgG
1903090	IgG

SHELF-LIFE : 1 year from date of shipment

PRESENTATION : Purified by chromatography by Protein G in PBS buffer, pH 7.4

IMMUNOGEN : Recombinant feline leukemia virus

CONCENTRATION : 1 - 2 mg/ml by OD at 280nm (extinction coefficient 1.4)

PRESERVATIVE : 0.1% sodium azide

QC ANALYSIS : By ELISA and SDS PAGE

APPLICATION : The antibodies are suitable for lateral flow and ELISA application using clone 1803080 as capture and clone 1903090 for detection.

STORAGE : $\leq -20^{\circ}\text{C}$.
Avoid repeated freezing and thawing

