



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: H01349M **Lot #:** 1C08117

Description: MAb to Hemoglobin
Monoclonal Antibody to Hemoglobin

Specificity: Hemoglobin

Host Animal: Mouse **Isotype:** IgG₁

Immunogen: Human Hemoglobin

Source: Ascites

Format: Purified, Liquid

Purification: > 95% pure (SDS-PAGE). Protein A Chromatography

Concentration: 5.37 mg/mL (OD280nm, E^{0.1%} = 1.4)

Buffer: 0.01 M Sodium Phosphate, 0.15 Sodium Chloride, pH 7.0.

Preservative: 0.02% Sodium Azide

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Recommended antibody pair for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
H01348M	H01349M

Storage: Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

CORRY THOMAS

Signature

22 March 2017

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY