

## Product Datasheet

### Heparin Binding Protein (HBP) Antibody HM722

<b>Product Name</b>	Heparin Binding Protein (HBP) Antibody
<b>Catalog Number</b>	HM722
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, CLIA
<b>Immunogen</b>	Purified Human Heparin Binding Protein
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	Heparin Binding Protein
<b>Matched Pairs</b>	HM721 (capture or detection)
<b>Preservatives</b>	0.05% Sodium Azide