

Product Datasheet

Recombinant IGF-1 antigen 610173

| | |
|--------------------------|------------------------------------------------------------------------------------------------|
| Product Name | Recombinant IGF-1 antigen |
| Catalog Number | 610173 |
| Description | Recombinant human IGF-1 antigen with a C-terminal histidine tag. Chromatographically purified. |
| Species/Host | E. coli |
| Alternative Names | Insulin-like Growth Factor 1, Human Insulin-like Growth Factor 1 |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid |
| Storage | -20°C |
| Buffer | 50 mM trisodium citrate, 0.9 % NaCl, pH 6.0 |
| Molecular Weight | Predicted molecular weight: 9 kDa |
| Shelf Life | Unspecified |