Dedication to the drug screening industry?

Drugs of Abuse catalog

Medix Biochemica



Introduction

Drugs of abuse is a broad term used for substances taken for reasons outside of the intended uses. Abuse of these substances often have a high potential of causing harm emotionally, physically and mentally on individuals, families and communities alike due to the mood- and mind-altering effects and dependencies that sometimes result. Testing for these substances (illegal or prescribed) is performed on individuals for several reasons. This may include treatment for substance use disorders, screening to ensure compliance in various programs and even assess criminal liability, among others. Various sample types, including blood, urine, saliva, sweat, and hair are used for a range of purposes and advantages. Tetrahydrocannabinol (THC), the major psychoactive ingredient in marijuana, and Ethyl glucuronide (EtG), an ethanol metabolite, are two common analytes tested by law enforcement or employment agencies to test intoxication or long-term use. Other causes for drugs of abuse testing

include medical screening, insurance purposes, and athletic testing among others. Specific analyte testing and broad screening panels are common. Medix Biochemica provides a wide selection of high-quality antigens and antibodies for drugs of abuse testing. Our optimized, industrial-scale manufacturing procedures, certified batch-to-batch consistency, and expert customer service have made Medix Biochemica an ideal raw material supply partner for assay and test kit manufacturers or other important customers in the in vitro diagnostic industry.

Products

6-MonoacetyImorphine (6-MAM)	12, 4
Amphetamine	4, 8
Benzodiazepine (Oxazepam)	8
Benzoylecgonine (BE)	4, 8
Buprenorphine (BUP)	8
Cannabinoid (THC)	8
Carisoprodol (Soma)	4, 9
Clonazepam	4
EDDP-BSA	9
Ethyl glucuronide (EtG)	5, 9
Fentanyl	5, 9
Hydrocodone	9
К2	10
Ketamine	5, 10
Kratom (Mitragynine)	10
Lysergic acid diethylamide (LSD)	5, 10
MDPV	11
Methadone	5
Methamphetamine	11
Methcathinone	11
Methylphenidate	6, 11
Morphine	6
Oxycodone	6
Phencyclidine (PCP)	6

, 4	Pinaca	12
, 8	Propoxyphene (PPX)	6
8	Tramadol	7
, 8	Tricyclic Antidepressant (TCA)	7, 12
8	UR-144	12
8	Zolpidem	7, 12

	Product	Concentration mg/mL *	Catalog number
Ethyl glucuronide (EtG)			
Ethyl glucuronide (EtG) is a metabolite of ethanol. It is used as a biomarker for alcohol consumption.	MONOCLONAL AN EtG R12401	TIBODY	140021 NEW
Fentanyl			
Fentanyl is a synthetic opioid with 50–100 times the potency of morphine. It is used to treat severe pain. Fentanyl may be illicitly produced and mixed with other illegal drugs increasing the abuse and overdose risks.	SYNTHETIC ANTIG Fentanyl-BSA Antigen	EN	170040 NEW
Ketamine			
Ketamine is a dissociative drug and anesthetic with hallucinogenic properties.	MONOCLONAL AN Ketamine Antibody	TIBODIES	100902 NEW
	SYNTHETIC ANTIGE	EN	
	Ketamine-BSA Antigen		170037 NEW
Lysergic acid diethylamide (I	LSD)		
Lysergic acid diethylamide (LSD) is a semi-synthetic hallucinogen.	MONOCLONAL AN LSD R12301	TIBODIES	TBD UPCOMIN
Methadone			
Methadone is a synthetic opioid agonist.	MONOCLONAL AN Methadone Antibody		100891 NEW

6-Monoacetylmorphine (6-MAM)

Heroin (diacetylmorphine) is metabolized to 6-monoacetylmorphine
and morphine. 6-monoacetylmorphine
is a unique metabolite of heroin and
its presence is definitive evidence of recent heroin use. Morphine is also
a metabolite of several legal drugs,
including codeine, and therefore is not a conclusive indication of heroin use.

SYNTHETIC ANTIGENS 6-Monoacetylmorphine (6-MAM)-BSA Antigen

Product

Concentration mg/mL *

170031 NEW

Catalog number

Amphetamine

Amphetamine is a central nervous	SYNTHETIC ANTIGENS	
system stimulant drug used to treat ADHD and narcolepsy.	Amphetamine-HRP Antigen	170035 NEW

Benzoylecgonine (BE)

Benzoylecgonine (BE) is the main	MONOCLONAL ANTIBODIES	
metabolite of cocaine.	Benzoylecgonine (BE) Antibody	100890 NEW

Carisoprodol (Soma)

Carisoprodol is a muscle relaxant used	SYNTHETIC ANTIGENS	
to treat skeletal muscle pain and injury.	Carisoprodol (Soma)-BSA Antigen	170046 NEW

Clonazepam

Clonazepam is a benzodiazepine that may
be used for the treatment of panic and
seizure disorders.

SYNTHETIC ANTIGEN Clonazepam-BSA Antigen

170064 NEW

Tramadol MONOCLONAL ANTIBODY Tramadol is a synthetic opioid drug used to treat pain. MONOCLONAL ANTIBODY Tramadol Antibody SYNTHETIC ANTIGEN Tramadol-BSA Tramadol-BSA Antigen Tramadol-BSA Difference Tramadol-BSA Antigen Tricyclic Antidepressant (TCA) Tricyclic antidepressants are a class of medications used primarily for the treatment of depression and other mood disorders. MONOCLONAL ANTIBODIES Tricyclic Anti-de-pressant (TCA) Tricyclic Anti-de-pressant (TCA) Antibody SYNTHETIC ANTIGEN		Product	Concentration mg/mL *
to treat pain. Tramadol Anti- body SYNTHETIC ANTIGEN Tramadol-BSA Antigen Tricyclic Antidepressant (TCA) Tricyclic antidepressants are a class of medications used primarily for the treatment of depression and other mood disorders. MONOCLONAL ANTIBODIES Tricyclic Anti-de- pressant (TCA) Antibody	Tramadol		
Tramadol-BSA Antigen Tricyclic Antidepressant (TCA) Tricyclic antidepressants are a class of medications used primarily for the treatment of depression and other mood disorders. MONOCLONAL ANTIBODIES		Tramadol Anti- body	
Tricyclic antidepressants are a class MONOCLONAL ANTIBODIES of medications used primarily for the Tricyclic Anti-de- treatment of depression and other pressant (TCA) mood disorders. Antibody			
Tricyclic Anti-		Antigen	

Zolpidem

Zolpidem is a central nervous system depressant drug causing drowsiness, sedation and decreased psychomotor function. Zolpidem is used to treat insomnia and may result in next-day impairment.

SYNTHETIC ANTIGEN Zolpidem-BSA Antigen

170041 NEW

Methylphenidate is a central nervous system stimulant drug used to treat Attention Deficient Hyperactivity Disorder (ADHD) and narcolepsy. Abuse risk increases when administered in excessive doses and via intravenous and intranasal routes. SYNTHETIC ANTIGENS Methylphenidate-BSA Antigen

Product

170039 NEW

Catalog number

Concentration mg/mL *

Morphine

Morphine is a non-synthetic opioid	MONOCLONAL ANTIBODIES	
drug used to treat pain.	Morphine Anti- body	100897 NEW
	SYNTHETIC ANTIGEN	
	Morphine-BSA Antigen	170036 NEW

Oxycodone

Oxycodone is a semi-synthetic opioid	MONOCLONAL ANTIBODIES	
drug used to treat pain.	Oxycodone Antibody	100899 NEW

Phencyclidine (PCP)

Phencyclidine is a dissociative	MONOCLONAL ANTIBODIES	
anesthetic that may cause hallucinations and violent behavior.	Phencyclidine (PCP) Antibody	100900 NEW

Propoxyphene (PPX)

Propoxyphene is an opioid used to treat	MONOCLONAL ANTIBODIES	
pain.	Propoxyphene (PPX) Antibody	100895 NEW
	Propoxyphene (PPX)-BSA Antigen	170032 NEW

Product	Source	Applications	Matched pairs	Catalog number
Amphetamine				
MONOCLONAL ANTIB	ODY			
Amphetamine Antibody	Mouse	ELISA, Lateral Flow	P91-99-71H (detection) or 170069 (capture-Espoo)	C03-99-72M-P
ANTIGEN Amphetamine-HRP Antigen		ELISA	C03-99-72M-P (detection)	P91-99-71H
Benzodiazepin	e (Oxazej	oam)		
POLYCLONAL ANTIBO	DY			
Benzodiazepine (Oxazepam) Antibody	Sheep	Lateral Flow	GN508 (capture)	GN100
ANTIGEN Benzodiazepine (Oxazepam)-BSA Antige	en	Lateral Flow	GN100 (detection)	GN508
Benzoylecgoni	ne (BE)			
ANTIGEN Benzoylecgonine (BE)- BSA Antigen		Lateral Flow	HM297 (detection)	LA004
Buprenorphine	e (BUP)			
ANTIGEN Buprenorphine (BUP)-BSA Antigen		ELISA, Lateral Flow	HM376 or HM377 (detection)	LA283
Cannabinoid (1	ГНС)			
POLYCLONAL ANTIBO	•			
Cannabinoid (THC) Antibody	Sheep	ELISA, Lateral Flow	LA020 (capture)	PP124IF

Product	Source	Applications	Matched pairs	Catalog number
Carisoprodol (S	Soma)			
POLYCLONAL ANTIBO	DY			
Carisoprodol (Soma) Antibody	Rabbit	ELISA	LA648 (capture) or 170046 (capture- Espoo)	P02-99-72R-IF
EDDP-BSA				
ANTIGEN EDDP-BSA Antigen		ELISA, Lateral Flow	HM381 (detection)	LA282
Ethyl glucuroni	ide (EtG)			
ANTIGEN Ethyl glucuronide (EtG)- BSA Antigen		ELISA, Lateral Flow	R12401/cat # 140021 (detection-Espoo)	LA189
Fentanyl				
POLYCLONAL ANTIBO	DY			
Fentanyl Antibody	Rabbit	ELISA	170040 (capture- Espoo) or P91-99-51H (detection)	P01-99-53R-IF
MONOCLONAL ANTIB	ODY			
Fentanyl Antibody	Mouse	Lateral Flow	LA008 (capture)	HM004
ANTIGEN				
Fentanyl-BSA Antigen Fentanyl-BSA Antigen		Lateral Flow ELISA, Lateral Flow	HM004 (detection) HM628 (detection)	LA008 LA395
Hydrocodone				
POLYCLONAL ANTIBO	DY			
Hydrocodone Antibody	Rabbit	ELISA, Lateral Flow	170054 (capture- Espoo)	JP101

Product	Source	Applications	Matched pairs	Catalog number
K2				
POLYCLONAL ANTIBO	DY			
K2 Antibody	Rabbit	ELISA, Lateral Flow	170082 (capture- Espoo)	JP099
MONOCLONAL ANTIB K2 Antibody	DDY Mouse	Lateral Flow	LA009 (capture)	HM005
Ketamine				
POLYCLONAL ANTIBO	DY			
Ketamine Antibody	Rabbit	ELISA	170037 (capture-Espoo)	PP211IF
MONOCLONAL ANTIB				
Ketamine Antibody Ketamine Antibody	Mouse Mouse	ELISA, Lateral Flow Lateral Flow	LA280 (capture) LA010 (capture)	HM390 HM006
	Mouse	Edicial How		110000
Ketamine-BSA Antigen		ELISA, Lateral Flow	HM390 or HM391 (detection)	LA280
Kratom (Mitrag	ynine)			
MONOCLONAL ANTIB	ODY			
Kratom (Mitragynine) Antibody	Mouse	ELISA, Lateral Flow		HM827
ANTIGEN Kratom (Mitragynine) Antigen		ELISA, Lateral Flow	HM1102 (detection)	LA653
Lysergic acid di	ethylamic	de (LSD)		
ANTIGEN Lysergic acid diethylamide (LSD)- BSA Antigen	A.	ELISA, Lateral Flow	HM1201 (detection)	LA783

Product	Source	Applications	Matched pairs	Catalog number
MDPV				
MONOCLONAL ANTI	BODY			
MDPV Antibody	Mouse	ELISA, Lateral Flow	LA011 (capture)	HM457
ANTIGEN MDPV-BSA Antigen		Lateral Flow	HM457 (detection)	LA011
Methampheta				
MONOCLONAL ANTI	BODY			
Methamphetamine Antibody	Mouse	Lateral Flow	P81-99-31A (capture) or P91-99-31H (detection)	C03-99-31M-P
Methamphetamine Antibody	Mouse	Lateral Flow	LA201 (capture)	HM165
Methcathinone Antibody	Mouse	ELISA		HM300
ANTIGEN				
Methamphetamine-BS Antigen	SA	Lateral Flow	HM165 (detection)	LA201

Methcathinone

MONOCLONAL ANTIE	ODY		
Methcathinone Antibody	Mouse	ELISA	HM300
ANTIGEN			
Methcathinone-BSA Antigen		Lateral Flow	LA267

Methylphenidate

POLYCLONAL ANTIBO	DY			
Methylphenidate Antibody	Rabbit	ELISA	170039 (capture- Espoo)	P03-99-02R-IF

Product	Source	Applications	Matched pairs	Catalog number
6-MonoacetyIn	orphine	(6-MAM)		
-				
POLYCLONAL ANTIBOI 6-Monoacetylmorphine (6-MAM) Antibody	Rabbit	ELISA, Lateral Flow	170031 (Espoo-capture)	JP092
MONOCLONAL ANTIBO	ODY			
6-MonoacetyImorphine (6-MAM) Antibody	Mouse	ELISA, Lateral Flow	170031 (Espoo-capture)	HM446
Pinaca				
MONOCLONAL ANTIB	ODY			
Pinaca Antibody	Mouse	ELISA, Lateral Flow	LA366 (capture)	HM548
Pinaca Antibody	Mouse	ELISA, Lateral Flow	LA384 (capture)	HM594
ANTIGEN			HM548, HM614, HM615	
		ELISA, Lateral Flow		LA366
Pinaca-BSA Antigen		ELISA, Edicitar I low	or HM616 (detection)	2,1000
Pinaca-BSA Antigen	pressant	ELISA, Lateral Flow	or HM616 (detection) HM594 (detection)	
•	-	ELISA, Lateral Flow		LA384
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant	DY	ELISA, Lateral Flow		LA384
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody	DY Sheep	ELISA, Lateral Flow		LA384
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-144	DY Sheep	ELISA, Lateral Flow		LA384 PP141IF
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-144 MONOCLONAL ANTIBO UR-144 Antibody ANTIGEN	Sheep	ELISA, Lateral Flow (TCA) ELISA ELISA	HM594 (detection)	LA384 PP141IF HM595
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-144 MONOCLONAL ANTIBO UR-144 Antibody	Sheep	ELISA, Lateral Flow (TCA) ELISA	HM594 (detection)	LA384 PP141IF HM595 LA385
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-144 MONOCLONAL ANTIBO UR-144 Antibody ANTIGEN	Sheep	ELISA, Lateral Flow (TCA) ELISA ELISA	HM594 (detection)	LA384 PP141IF HM595
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-144 MONOCLONAL ANTIBO UR-144-Antibody ANTIGEN UR-144-BSA Antigen Zopiclone ANTIGEN	Sheep	ELISA, Lateral Flow (TCA) ELISA ELISA, Lateral Flow ELISA, Lateral Flow	HM594 (detection)	LA384 PP141IF HM595 LA385
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-1444 MONOCLONAL ANTIBO UR-144-BSA Antigen Zopiclone ANTIGEN Zopiclone	DY Sheep DDY Mouse	ELISA, Lateral Flow (TCA) ELISA ELISA	HM594 (detection)	LA384 PP141IF HM595
Pinaca-BSA Antigen Tricyclic Antide POLYCLONAL ANTIBOI Tricyclic Antidepressant (TCA) Antibody UR-144 MONOCLONAL ANTIBO UR-144-Antibody ANTIGEN UR-144-BSA Antigen Zopiclone ANTIGEN	DY Sheep DDY Mouse	ELISA, Lateral Flow (TCA) ELISA ELISA, Lateral Flow ELISA, Lateral Flow	HM594 (detection)	LA384 PP141IF HM595 LA385

Medix Biochemica Group

Medix Biochemica is a global, market-leading supplier of critical raw materials to the *in vitro* diagnostics (IVD) industry. Our mission is to enable IVD customers to deliver accurate results to billions of patients worldwide. By working in partnership with Medix Biochemica, IVD manufacturers can reduce their time to market through our industry-leading expertise which increases efficiencies and builds quality into their tests. We achieve this through our deep scientific know-how combined with the broadest and most relevant raw material portfolio and our relentless focus on quality.

Medix Biochemica Group employs more than 260 professionals worldwide. With headquarters in Finland, we serve our worldwide customers through our local teams in the US, Europe, and Asia. www.medixbiochemica.com



The first-choice partner

- Focusing on IVD manufacturers' needs
- Capability to support customers across the globe
- Focus on new product development based on customer needs
- Efficient and scalable manufacturing capabilities



Widest portfolio

- Monoclonal antibodies
- Antigens for immunoassay materials
- Molecular diagnostics reagents
- Enzymes for clinical chemistry controls
- Biochemicals and biologicals



Highest quality

- ISO 13485 & FDA QSR compliance
- Batch to batch consistency
- Scalable production capability from milligrams to kilograms

Close to our global customers

We serve our global customers from our sites in Finland, the USA, France, Germany and China.



Our values



We care

We put the customer first; we care for the people and the world around us.



We dare

We are ambitious and innovative; we take risks and dare to win with our customers.

We deliver

We commit to our customers and keep our promises.

Get in touch

Orders: orders@medixbiochemica.com

Sales inquiries: sales@medixbiochemica.com

Molecular Diagnostics: MDx@medixbiochemica.com

Technical support: techsupport@medixbiochemica.com

General inquiries: medix@medixbiochemica.com

Medix Biochemica

Ordering and product information **medixbiochemica.com**

Copyright © 07/2022 Medix Biochemica. All rights reserved.

Lee Biosolutions, EastCoast Bio, Diaclone, myPOLS Biotech, and Bioresource Technology are part of Medix Biochemica Group. Medix Biochemica reserves the right to make changes and improvements to any of the products described in this document without prior notice.