

Product Datasheet

Pregnancy-Associated Plasma Protein-A (PAPP-A) 496-90

Product Name	Pregnancy-Associated Plasma Protein-A (PAPP-A)
Catalog Number	496-90
Alternative Names	PAPP-A, Pappalysin-1
Brand	Lee Biosolutions
Form/Appearance	Frozen Liquid in Phosphate Buffered Saline, pH 7.4
Concentration	≥ 1.0 mg/mL (Coomasie)
Storage	-20 °C
Purity	≥ 85% (SDS-PAGE)
Source	Human Serum
Associated Products	PAPP-A Monoclonal Antibody HM1014 and HM1015
Molecular Weight	~200,000
Shelf Life	2 years