

Product Datasheet

Anti-h ProGRP 13001 SPTN-5 100989

Product Name	Anti-h ProGRP 13001 SPTN-5
Catalog Number	100989
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Alternative Names	progastrin-releasing peptide
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	aa 45-46
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.1-6.7
Specificity	Antibody recognizes human progastrin-releasing peptide
Matched Pairs	Anti-h ProGRP 13002, HM784 and HM785
Shelf Life	Unspecified